

5 November 2013 EMA/675665/2013 Information and Communications Technology

EudraCT v9.0 Release Notes

EudraCT Main version: 8.4.9.1

EudraCT Results version: 9.0.1.1

Goal: Patch release

Release Date: 09/11/2013

Release Notes for EudraCT v9.0.1 Public

The functionality offered by EudraCT v9.0.1 is the same as v9.0.0.

The following problems have been fixed since the last release:

- In the Results module, warning messages for serious and non-serious adverse events are displayed when values for subjects exposed are equal to the total number of subjects exposed for a reporting group.
- After a user has pressed the 'confirm' button to update a finalised result, the application sometimes takes a while to proceed to the next step. If the user presses the 'confirm' button again while waiting for the application to process the confirm request, the following error message is displayed: "An unexpected application error has occurred".
- The content of the eudract.ema.europa.eu website has been updated
- The Help has been updated with additional and improved content.

Known issues in release 9.0.1

- When the user of the Results module attaches several summary files to the result, in certain cases an error message is displayed. If this happens, the user is advised to go back and to repeat the attachment procedure.
- It may take several seconds for EudraCT to switch from the Results module to the Clinical Trial Application module. This delay is observable for the first time only, subsequent switches are quicker.
- Age continuous list of units is not standardised with EUTCT.
- There are a few tool tips missing or misplaced within the Results module.

7 Westferry Circus • Canary Wharf • London E14 4HB • United Kingdom **Telephone** +44 (0)20 7418 8400 **Facsimile** +44 (0)20 7418 8669 **E-mail** info@ema.europa.eu **Website** www.ema.europa.eu



An agency of the European Union

© European Medicines Agency, 2013. Reproduction is authorised provided the source is acknowledged.

• Some of the graphical elements are misplaced in certain screen resolutions in the CTA module due to the refreshed appearance in v9.0.0.